

Bis(4-cyanophenyl)methanol


Catalog No: #S0534

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bis(4-cyanophenyl)methanol
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	234.25
Formulation	C15H10N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(C1=CC=C(C=C1)C#N)C1=CC=C(C=C1)C#N
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Precht JC, et al. Drug Metab Dispos. 2013 Nov;41(11):1906-13.

Note: This product is for in vitro research use only