

## Cyclopropanecarbonitrile

Catalog No: #S0536


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cyclopropanecarbonitrile
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	67.09
Formulation	C4H5N
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular SystemSMILES:N#CC1CC1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Garca Ruano JL, et al. Org Lett. 2004 Dec 23;6(26):4945-8.

Note: This product is for in vitro research use only