

Desbenzyl Donepezil Hydrochloride

Catalog No: #S0537


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Desbenzyl Donepezil Hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	325.83
Formulation	C17H24ClNO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: Cl.COC1=C(OC)C=C2C(=O)C(CC3CCNCC3)CC2=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited solubility
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ferrarese A, et al. Xenobiotica. 1993 Sep;23(9):1007-16.

Note: This product is for in vitro research use only