

Dibenzosuberol

Catalog No: #S0540


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dibenzosuberol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	210.28
Formulation	C15H14O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:OC1C2=CC=CC=C2CCC2=CC=CC=C12Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Murray MA, Babe LM. Antiviral Res. 1999 Dec 15;44(2):123-31.

Note: This product is for in vitro research use only