

Fluoronaphthalene

Catalog No: #S0541


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fluoronaphthalene
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	146.16
Formulation	C10H7F
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:FC1=CC=CC2=CC=CC=C12Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Cerniglia CE, et al. Appl Environ Microbiol. 1984 Aug;48(2):294-300.

Note: This product is for in vitro research use only