

## Homoisovanillic acid


Catalog No: #S0542

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Homoisovanillic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	182.17
Formulation	C <sub>9</sub> H <sub>10</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:COC1=C(O)C=C(CC(O)=O)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Mathieu P, et al. J Neurochem. 1972 Jan;19(1):81-6.

Note: This product is for in vitro research use only