

Dichlorodiphenylmethane

Catalog No: #S0544


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dichlorodiphenylmethane
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	237.12
Formulation	C13H10Cl2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:ClC(Cl)(C1=CC=CC=C1)C1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Reid TM, et al. Toxicol Lett. 1998 May; 95(3):205-10.

Note: This product is for in vitro research use only