

## Methyl 5-methylnicotinate

Catalog No: #S0548


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Methyl 5-methylnicotinate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	151.16
Formulation	C8H9NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: COC(=O)C1=CC(C)=CN=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Limited  
solubility  
Appearance: White Needle Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Katzman M, et al. Neuropsychopharmacology. 2003 Oct; 28(10): 1846-51.

Note: This product is for in vitro research use only