

4'-tert-Butylacetophenone


Catalog No: #S0549

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4'-tert-Butylacetophenone
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	176.25
Formulation	C12H16O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC(=O)C1=CC=C(C=C1)C(C)(C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Ethanol, insoluble in WaterAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only