

Isonicotinamide

Catalog No: #S0552


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Isonicotinamide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	122.12
Formulation	C6H6N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: NC(=O)C1=CC=NC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Water: 191.7 g/L (37 °C)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Fukaya M, et al. Biochem Biophys Res Commun. 2013 Dec 6;442(1-2):92-8.

Note: This product is for in vitro research use only