

4-Hydroxyisophthalic acid


Catalog No: #S0553

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Hydroxyisophthalic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	182.13
Formulation	C8H6O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(=O)C1=CC(C(O)=O)=C(O)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Piscopo E, et al. Boll Soc Ital Biol Sper. 1984 Jun 30;60(6):1169-75.

Note: This product is for in vitro research use only