

Dibutylurea

Catalog No: #S0554


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dibutylurea
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	172.27
Formulation	C ₉ H ₂₀ N ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCCCNC(=O)NCCCC
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Off-white Color Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Bischoff M, et al. Biodegradation. 2005 Jun;16(3):265-73.

Note: This product is for in vitro research use only