

4-Diethylamino-2-butyne-1-ol


Catalog No: #S0555

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Diethylamino-2-butyne-1-ol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	141.21
Formulation	C ₈ H ₁₅ NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCN(CC)CC#CCO
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only