

3-(Methylamino)toluene

Catalog No: #S0557


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	3-(Methylamino)toluene
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	121.18
Formulation	C8H11N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CNC1=CC=CC(C)=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Patil ST, et al. Fresenius J Anal Chem. 2000 May;367(1):91-3.

Note: This product is for in vitro research use only