

Hydroferulic acid

Catalog No: #S0558


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hydroferulic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	196.2
Formulation	C10H12O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=C(O)C=CC(CCC(O)=O)=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Garcia-Conesa MT, et al. Redox Rep. 1997 Oct-Dec;3(5-6):319-23.

Note: This product is for in vitro research use only