

## Dibenzosuberone

Catalog No: #S0559


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dibenzosuberone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	208.26
Formulation	C15H12O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:O=C1C2=CC=CC=C2CCC2=CC=CC=C12Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Merkas S, et al. Molecules. 2005 Dec 31;10(12):1429-37.

Note: This product is for in vitro research use only