

## 3,3-Diphenylpropylamine


Catalog No: #S0560

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	3,3-Diphenylpropylamine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	211.3
Formulation	C15H17N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular System  
SMILES: NCCC(C1=CC=CC=C1)C1=CC=CC=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting  
pt: Solubility: Limited solubility  
Appearance: Solid-Liquid Mixture  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. KOCHSIEK K, et al. *Arzneimittelforschung*. 1960 Aug;10:576-83.

Note: This product is for in vitro research use only