

Protocatechuic acid

Catalog No: #S0562


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protocatechuic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	154.12
Formulation	C7H6O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:OC(=O)C1=CC(O)=C(O)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lv Y, et al. Biotechnol Lett. 2014 Nov; 36(11):2303-10.

Note: This product is for in vitro research use only