

Homoveratronitrile

Catalog No: #S0563


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Homoveratronitrile
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	177.2
Formulation	C10H11NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=C(OC)C=C(CC#N)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Grodetskaia NS, et al. Gig Tr Prof Zabol. 1982 Feb;(2):55.

Note: This product is for in vitro research use only