

1-Indanone

Catalog No: #S0564


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	1-Indanone
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	132.16
Formulation	C ₉ H ₈ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: O=C1CCC2=CC=CC=C12
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Water: 6.5 g/L at 20 °C
Appearance: White Sheet Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Fill TP, et al. J Microbiol Biotechnol. 2012 Jun;22(6):832-7.

Note: This product is for in vitro research use only