

Bis(2-chloroethyl)amine hydrochloride

Catalog No: #S0565


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bis(2-chloroethyl)amine hydrochloride
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	178.49
Formulation	C4H9CL2NHCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:Cl.CCCNCCCIPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:Off-white
Color PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can
be stored below -20 °C for several month.

References

1. HUNT CC, PHILIPS FS. J Pharmacol Exp Ther. 1949 Feb;95(2):131-44.

Note: This product is for in vitro research use only