

2-Methylbenzhydrol

Catalog No: #S0566

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2-Methylbenzhydrol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	198.26
Formulation	C14H14O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC1=CC=CC=C1C(O)C1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Bose KS, Sarma RH. Biochem Biophys Res Commun. 1975 Oct 27;66(4):1173-9.

Note: This product is for in vitro research use only