

Homoveratrumic acid

Catalog No: #S0568


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Homoveratrumic acid |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 196.2 |
| Formulation | C10H12O4 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:Nervous systemSMILES:COC1=C(OC)C=C(CC(O)=O)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White Needle CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Dana M, et al. Molecules. 2011 May 5;16(5):3826-44.

Note: This product is for in vitro research use only