

## Pyrazinoic acid

Catalog No: #S0569


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pyrazinoic acid
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	124.1
Formulation	C5H4N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:OC(=O)C1=NC=CN=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in hot WaterAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Rada B, et al. Acta Virol. 1971 Jul;15(4):326-8.

Note: This product is for in vitro research use only