

3,5-Dihydroxybenzoic acid

Catalog No: #S0570


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | 3,5-Dihydroxybenzoic acid |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Lipid |
| Calculated MW | 154.12 |
| Formulation | C7H6O4 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
SMILES: OC(=O)C1=CC(O)=CC(O)=C1
Pathways: Metabolism
Receptor: lipolysis
Boiling pt: Melting pt: Solubility: Water: 84 g/L (20 C)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Liu C, et al. J Pharmacol Exp Ther. 2012 Jun;341(3):794-801.

Note: This product is for in vitro research use only