

2-Methylbenzene-1-sulfonamide

Catalog No: #S0571


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2-Methylbenzene-1-sulfonamide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	171.22
Formulation	C7H9NO2S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC1=CC=CC=C1S(N)(=O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kacprzak JL. J Assoc Off Anal Chem. 1978 Nov;61(6):1528-32.

Note: This product is for in vitro research use only