

2-Methoxynaphthalene

Catalog No: #S0573


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2-Methoxynaphthalene
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	158.2
Formulation	C11H10O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=CC2=CC=CC=C2C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Mohapatra H, Umapathy S. J Phys Chem A. 2010 Dec 2;114(47):12447-51.

Note: This product is for in vitro research use only