

1-Methylimidazole


Catalog No: #S0575

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	1-Methylimidazole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	82.1
Formulation	C4H6N2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CN1C=CN=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:soluble in

WaterAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Shafizadeh N, et al. J Chem Phys. 2014 Nov 7;141(17):174310.

Note: This product is for in vitro research use only