

Diethylaminoethoxyethanol


Catalog No: #S0576

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Diethylaminoethoxyethanol
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	161.24
Formulation	C8H19NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCN(CC)CCOCCO
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. CUTLER A, et al. Endocrinology. 1961 Sep;69:473-82.

Note: This product is for in vitro research use only