

## Erythrosin B

Catalog No: #S0578


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Erythrosin B
Brief Description	Inhibitors
Purification	98.36%
Target Name	Others
Calculated MW	835.89
Formulation	C20H8I4O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: OC1=C(I)C=C2C(OC3=C(I)C(O)=C(I)C=C3C22OC(=O)C3=CC=CC=C23)=C1I  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: Red Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Horikoshi S, et al. Clin Chim Acta. 2012 Jul 11;413(13-14):1087-91.

Note: This product is for in vitro research use only