

## 4-Pentynoic acid


Catalog No: #S0580

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	4-Pentynoic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	98.1
Formulation	C <sub>5</sub> H <sub>6</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: OC(=O)CCC#C Pathways: Others Receptor: Others Boiling pt: 54-57 Melting pt: Solubility: Soluble in DMSO Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only