

Erythritol

Catalog No: #S0583


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Erythritol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	122.12
Formulation	C4H10O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:C([C@H]([C@H](CO)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water: 24 mg/mL (196.52 mM); DMSO: 24 mg/mL (196.52 mM)Appearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Munro IC, et al. Food Chem Toxicol, 1998, 36(12), 1139-1174.

Note: This product is for in vitro research use only