

D-(+)-Maltose monohydrate

Catalog No: #S0584L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	D-(+)-Maltose monohydrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	360.31
Formulation	C ₁₂ H ₂₂ O ₁₁ H ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C([C@@H]1[C@H]([C@@H]([C@H]([C@H](O1)O)[C@H]([C@@H](CO)O)[C@@H]([C@H](C=O)O)O)O)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: 119-121°C Solubility: Soluble in Water Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sugawara K, et al. Anal Chim Acta. 2014 Mar 3;814:55-62.

Note: This product is for in vitro research use only