

## L(+)-Arabinose

Catalog No: #S0587


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	L(+)-Arabinose
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	150.13
Formulation	C5H10O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
SMILES: C1[C@@H]([C@@H]([C@H](C(O1)O)O)O)O  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Soluble in Water  
Appearance: White Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Rodrigues Borges M, Balaban Rde C. J Biotechnol. 2014 Dec 20;192 Pt A:42-9.

Note: This product is for in vitro research use only