

## Phenylhydantoin


Catalog No: #S0596

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phenylhydantoin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	176.17
Formulation	C9H8N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:O=C1NC(C(=O)N1)C1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Serra R, et al. J Inflamm (Lond). 2010 Sep 15;7:48.

Note: This product is for in vitro research use only