

N-Methyl-N-(trimethylsilyl)trifluoroacetamide


Catalog No: #S0599

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-Methyl-N-(trimethylsilyl)trifluoroacetamide
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	199.25
Formulation	C ₆ H ₁₂ F ₃ NOSi
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CN(C(=O)C(F)(F)F)[Si](C)(C)C
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhou YQ, et al. J Environ Sci (China). 2007;19(7):879-84.

Note: This product is for in vitro research use only