

## Phthalic Acid Monobenzyl Ester


Catalog No: #S0601

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phthalic Acid Monobenzyl Ester
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	256.25
Formulation	C15H12O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: OC(=O)C1=CC=CC=C1C(=O)OCC1=CC=CC=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting  
pt: Solubility: Limited solubility  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Serrano SE, et al. Int J Environ Res Public Health. 2014 Jun 12;11(6):6193-215.

Note: This product is for in vitro research use only