

## Phthalic Acid Monooctyl Ester


Catalog No: #S0603

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phthalic Acid Monooctyl Ester
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	278.34
Formulation	C16H22O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:CCCCC(CC)COC(=O)C1=CC=CC=C1C(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Ellero-Simatos S, et al. J Proteome Res. 2011 Dec 2;10(12):5493-502.

Note: This product is for in vitro research use only