

Dihydroartemisinin

Catalog No: #S0607


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dihydroartemisinin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	284.35
Formulation	C15H24O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C[C@@H]1CC[C@H]2[C@H]([C@H](O[C@H]3[C@@]24[C@H]1CCC(O3)(OO4)C)O)C Pathways: Others Receptor: Others
 oling pt: Melting pt: Solubility: Ethanol: 11 mg/mL warmed (38.68 mM); DMSO: 56 mg/mL warmed (196.94 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Yang D, et al. Int J Clin Exp Pathol. 2015 Feb 1;8(2):1270-81.

Note: This product is for in vitro research use only