

## 4-Hydroxy-3-methoxymandelic acid


Catalog No: #S0611

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	4-Hydroxy-3-methoxymandelic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	198.17
Formulation	C9H10O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:COC1=C(O)C=CC(=C1)C(O)C(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Kolentinis MK, et al. Eur J Pharmacol. 2013 Mar 5;703(1-3):46-52.

Note: This product is for in vitro research use only