

Beta-Hydroxyisovaleric acid

Catalog No: #S0613


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-Hydroxyisovaleric acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	118.13
Formulation	C5H10O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC(C)(O)CC(=O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited
solubilityAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a
while.Stock solution can be stored below -20 °C for several month.

References

1. Sabourin PJ, Bieber LL. Arch Biochem Biophys. 1981 Jan;206(1):132-44.

Note: This product is for in vitro research use only