

1-(3-Chlorophenyl)piperazine hydrochloride

Catalog No: #S0615


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	1-(3-Chlorophenyl)piperazine hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	233.14
Formulation	C10H14Cl2N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: Cl.C1C=CC(=C1)N1CCNCC1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Gray Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Arbo MD, et al. Drug Alcohol Depend. 2012 May 1;122(3):174-85.

Note: This product is for in vitro research use only