

Indole-3-Glyoxylyl Chloride

Catalog No: #S0618


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|-------------------------------------------------|
| Product Name | Indole-3-Glyoxylyl Chloride |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Autophagy activator |
| Calculated MW | 207.62 |
| Formulation | C10H6ClNO2 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Immunology
SMILES: ClC(=O)C(=O)C1=CNC2=CC=CC=C12
Pathways: Autophagy Receptor; Lysozyme
Boiling pt: Melting pt: Solubility: Limited solubility
Appearance: Brown Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Nakazono M, et al. Anal Sci. 2003 Jan;19(1):123-7.

Note: This product is for in vitro research use only