

## Diphenyl oxide


Catalog No: #S0621

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Diphenyl oxide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	170.21
Formulation	C12H10O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: O(C1=CC=CC=C1)C1=CC=CC=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Water: 18 mg/l  
Appearance: White Needle Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Garcia-Domenech R, et al. Mol Divers. 2015 May;19(2):357-66.

Note: This product is for in vitro research use only