

## Diphenyl oxide

Catalog No: #S0621


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Diphenyl oxide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	170.21
Formulation	C12H10O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:O(C1=CC=CC=C1)C1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water: 18 mg/l  
Appearance:White Needle CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Garcia-Domenech R, et al. Mol Divers. 2015 May;19(2):357-66.

Note: This product is for in vitro research use only