

4-Hydroxybenzophenone

Catalog No: #S0624


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4-Hydroxybenzophenone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	198.22
Formulation	C13H10O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: OC1=CC=C(C=C1)C(=O)C2=CC=CC=C2 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited solubility Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Buck Louis GM, et al. Am J Epidemiol. 2014 Dec 15;180(12):1168-75.

Note: This product is for in vitro research use only