

Uric Acid

Catalog No: #S0626


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Uric Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Phosphorylase agonist
Calculated MW	168.11
Formulation	C ₅ H ₄ N ₄ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: O=C1NC2=C(N1)C(=O)NC(=O)N2
Pathways: Others
Receptor: Glycogen phosphorylase
Boiling pt: Melting
pt: Solubility: Water: 60 mg/L (at 20 °C)
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Livingston JR, et al. J Obstet Gynaecol Can. 2014 Oct;36(10):870-7.

Note: This product is for in vitro research use only