

beta-Alanine


Catalog No: #S0627

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	beta-Alanine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Amino Acids and Derivatives
Calculated MW	89.09
Formulation	C3H7NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:NCCC(O)=OPathways:OthersReceptor:Amino Acids and DerivativesBoiling pt:Melting pt:Solubility:soluble in WaterAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Hill CA, et al. Amino Acids. 2007 Feb;32(2):225-33.

Note: This product is for in vitro research use only