

## Atropic acid

Catalog No: #S0629


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Atropic acid
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	148.16
Formulation	C <sub>9</sub> H <sub>8</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:OC(=O)C(=C)C1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. GOTTLIEB S. J Am Chem Soc. 1948 Jan;70(1):423.

Note: This product is for in vitro research use only