

Ethyleneurea


Catalog No: #S0630

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ethyleneurea
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	86.09
Formulation	C3H6N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: O=C1NCCN1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: soluble in Water
Appearance: White
Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Nilsson BM, et al. J Med Chem. 1992 Aug 21;35(17):3270-9.

Note: This product is for in vitro research use only